

For Research Use Only

Nicalin Polyclonal antibody

Catalog Number: 18444-1-AP



Basic Information

Catalog Number:

18444-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13318

GenBank Accession Number:

BC010064

GeneID (NCBI):

56926

UNIPROT ID:

Q969V3

Full Name:

nicalin homolog (zebrafish)

Calculated MW:

63 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Nicalin, also named as NCLN, is distantly related to Nicastrin, a component of the Alzheimer's disease-associated γ -secretase, and forms a complex with Nomo. It also forms a complex with NOMO2 and TMEM147.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

