

For Research Use Only

# IGFBP6 Polyclonal antibody

Catalog Number: 18445-1-AP **1 Publications**



## Basic Information

Catalog Number: 18445-1-AP	GenBank Accession Number: BC011708	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 3489	
Source: Rabbit	Full Name: insulin-like growth factor binding protein 6	
Isotype: IgG	Calculated MW: 25 kDa	
Immunogen Catalog Number: AG13322		

## Applications

Tested Applications:  
ELISA

Cited Applications:  
WB

Species Specificity:  
human, mouse, rat

Cited Species:  
human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Alexandru Tatomir	33569054	Front Immunol	WB

## Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data