

For Research Use Only

# HOXD8 Polyclonal antibody

Catalog Number: 18588-1-AP



## Basic Information

### Catalog Number:

18588-1-AP

### Size:

150ul , Concentration: 300 ug/ml by Nanodrop and 260 ug/ml by Bradford method using BSA as the standard;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG13245

### GenBank Accession Number:

BC038709

### GeneID (NCBI):

3234

### UNIPROT ID:

P13378

### Full Name:

homeobox D8

### Calculated MW:

32 kDa

### Observed MW:

32 kDa

### Purification Method:

Antigen affinity purification

### Recommended Dilutions:

WB 1:200-1:1000

## Applications

### Tested Applications:

WB, ELISA

### Species Specificity:

human, mouse, rat

### Positive Controls:

WB : A375 cells, HEK-293 cells, HeLa cells

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



A375 cells were subjected to SDS PAGE followed by western blot with 18588-1-AP (HOXD8 Antibody) at dilution of 1:300 incubated at 4 degree celsius over night.