

For Research Use Only

EZH2 Polyclonal antibody

Catalog Number: 18691-1-AP **1 Publications**



Basic Information

Catalog Number: 18691-1-AP	GenBank Accession Number: BC010858	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 187 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 2146	
Source: Rabbit	UNIPROT ID: Q15910	
Isotype: IgG	Full Name: enhancer of zeste homolog 2 (Drosophila)	
Immunogen Catalog Number: AG6556	Calculated MW: 751 aa, 86 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse, rat

Notable Publications

Author	Pubmed ID	Journal	Application
Yuelong Jiang	36355485	Pharmaceuticals (Basel)	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data