

For Research Use Only

GABRA6 Polyclonal antibody

Catalog Number: 18728-1-AP



Basic Information

Catalog Number: 18728-1-AP	GenBank Accession Number: BC096241	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 2559	
Source: Rabbit	Full Name: gamma-aminobutyric acid (GABA) A receptor, alpha 6	
Isotype: IgG	Calculated MW: 453 aa, 51 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

GABRA6 belongs to the ligand-gated ionic channel (TC 1.A.9) family. GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data