

For Research Use Only

GABRR2 Polyclonal antibody

Catalog Number: 18730-1-AP



Basic Information

Catalog Number:

18730-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC130352

GeneID (NCBI):

2570

UNIPROT ID:

P28476

Full Name:

gamma-aminobutyric acid (GABA)
receptor, rho 2

Calculated MW:

57 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

Background Information

GABRR2 belongs to the ligand-gated ionic channel (TC 1.A.9) family. GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel. GABRR2 plays a role in retinal neurotransmission.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

