

For Research Use Only

HOXB3 Polyclonal antibody

Catalog Number: 18733-1-AP



Basic Information

Catalog Number:

18733-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002146

GeneID (NCBI):

3213

UNIPROT ID:

P14651

Full Name:

homeobox B3

Calculated MW:

44 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

HOXB3, also named as HOX2G, belongs to the Antp homeobox family. HOXB3 is sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. HOXB3 plays a role in the proliferation and differentiation of both early myeloid and lymphoid developmental pathways.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

