

For Research Use Only

HOXB3 Polyclonal antibody

Catalog Number: 18733-1-AP



Basic Information

Catalog Number: 18733-1-AP	GenBank Accession Number: NM_002146	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 3213	
Source: Rabbit	Full Name: homeobox B3	
Isotype: IgG	Calculated MW: 44 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

HOXB3, also named as HOX2G, belongs to the Antp homeobox family. HOXB3 is sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. HOXB3 plays a role in the proliferation and differentiation of both early myeloid and lymphoid developmental pathways.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data