

For Research Use Only

WNT8A Polyclonal antibody

Catalog Number: 18740-1-AP



Basic Information

Catalog Number:

18740-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC156844

GeneID (NCBI):

7478

UNIPROT ID:

Q9H1J5

Full Name:

wingless-type MMTV integration site family, member 8A

Calculated MW:

39 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

Background Information

WNT8A, also named as WNT8D, belongs to the Wnt family of secreted signaling proteins. It is a ligand for members of the frizzled family of seven transmembrane receptors. These have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. WNT8A may play an important role in the development and differentiation of certain forebrain structures, notably the hippocampus.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

