

For Research Use Only

# ABCA8 Polyclonal antibody

Catalog Number: 18741-1-AP



## Basic Information

<b>Catalog Number:</b> 18741-1-AP	<b>GenBank Accession Number:</b> bc130280	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul	<b>GeneID (NCBI):</b> 10351	
<b>Source:</b> Rabbit	<b>Full Name:</b> ATP-binding cassette, sub-family A (ABC1), member 8	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 1581 aa, 179 kDa	

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human

## Background Information

ABCA8 belongs to the ATP binding cassette (ABC) transporter family and it is a member of the ABC1 subfamily. ABC proteins transport various molecules across extra and intracellular membranes. The function of ABCA8 has not yet been determined.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data