

For Research Use Only

# ADRA1D Polyclonal antibody, PBS Only

Catalog Number: 18823-1-PBS



## Basic Information

Catalog Number:

18823-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_000678

GeneID (NCBI):

146

UNIPROT ID:

P25100

Full Name:

adrenergic, alpha-1D-, receptor

Calculated MW:

60 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

ADRA1D, also known as ADRA1A, belongs to the G-protein-coupled receptor 1 family. ADRA1D is involved in multicellular organism aging, neuron-glia cell signaling, and regulation of blood pressure.

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 18823-1-AP (ADRA1D antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 18823-1-PBS in a different storage buffer formulation.