

For Research Use Only

# HMCN1 Polyclonal antibody, PBS Only

Catalog Number:18837-1-PBS



## Basic Information

**Catalog Number:**

18837-1-PBS

**Size:**

100ug , Concentration: 1 mg/ml by Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**GenBank Accession Number:**

NM\_031935

**GeneID (NCBI):**

83872

**UNIPROT ID:**

Q96RW7

**Full Name:**

hemicentin 1

**Calculated MW:**

613 kDa, 600 kDa, 60 kDa

**Observed MW:**

60 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

WB, IHC, IF/ICC, Indirect ELISA

**Species Specificity:**

human, mouse, rat

## Background Information

Hemicentin-1 (HMCN1), also known as fibulin-6 (FIBL-6), is an extracellular protein of the immunoglobulin superfamily. HMCN1 consists of a von Willebrand A domain, tandem immunoglobulin modules, tandem EGF modules, and a fibulin carboxy-terminal module (PMID: 24912920). Mutations of HMCN1 gene may be associated with age-related macular degeneration. Three alternative transcripts encoding different HMCN1 isoforms have been described. This antibody raised against a synthesized peptide corresponding to 960-976 amino acids of human HMCN1 can recognize all three isoforms of HMCN1.

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human colon using 18837-1-AP (HMCN1 antibody) at dilution of 1:50 (under 10x lens). This data was developed using the same antibody clone with 18837-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon using 18837-1-AP (HMCN1 antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 18837-1-PBS in a different storage buffer formulation.



mouse skin tissue were subjected to SDS PAGE followed by western blot with 18837-1-AP (HMCN1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 18837-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed hTERT-RPE1 cells using HMCN1 antibody (18837-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 18837-1-PBS in a different storage buffer formulation.