

For Research Use Only

# UBN2 Polyclonal antibody, PBS Only

Catalog Number: 18844-1-PBS



## Basic Information

Catalog Number:

18844-1-PBS

Size:

100ug, Concentration: 1 mg/ml by  
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_173569

GeneID (NCBI):

254048

UNIPROT ID:

Q6ZU65

Full Name:

ubiquitin 2

Calculated MW:

146 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human testis using 18844-1-AP (UBN2 antibody) at dilution of 1:50 (under 10x lens). This data was developed using the same antibody clone with 18844-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human testis using 18844-1-AP (UBN2 antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 18844-1-PBS in a different storage buffer formulation.