

For Research Use Only

MAP3K15 Polyclonal antibody

Catalog Number: 18849-1-AP



Basic Information

Catalog Number: 18849-1-AP	GenBank Accession Number: NM_001001671	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 389840	
Source: Rabbit	Full Name: mitogen-activated protein kinase kinase kinase 15	
Isotype: IgG	Calculated MW: 147 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

MAP3K15 belongs to the protein kinase superfamily; STE Ser/Thr protein kinase family and MAP kinase kinase kinase subfamily. It is a component of a protein kinase signal transduction cascade. It contains 1 protein kinase domain.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data