

For Research Use Only

# GJE1 Polyclonal antibody

Catalog Number: 18851-1-AP



## Basic Information

Catalog Number: 18851-1-AP	GenBank Accession Number: AL033522	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 100126572	
Source: Rabbit	Full Name: gap junction protein, epsilon 1, 23kDa	
Isotype: IgG	Calculated MW: 24 kDa	

## Applications

Tested Applications:  
ELISA  
Species Specificity:  
human, mouse, rat

## Background Information

GJE1, also named as Cx23, belongs to the connexin family and Beta-type (group I) subfamily. It is a connexon which is composed of a hexamer of connexins

## Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data