

For Research Use Only

AKR1C4 Polyclonal antibody, PBS Only

Catalog Number: 18852-1-PBS



Basic Information

Catalog Number:
18852-1-PBS

Size:
100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:
Rabbit

Isotype:
IgG

Immunogen Catalog Number:
AG13444

GenBank Accession Number:
BC020744

GeneID (NCBI):
1109

UNIPROT ID:
P17516

Full Name:
aldo-keto reductase family 1,
member C4 (chlordecone reductase;
3-alpha hydroxysteroid
dehydrogenase, type I; dihydrodiol
dehydrogenase 4)

Calculated MW:
323 aa, 37 kDa

Observed MW:
37 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 18852-1-AP (AKR1C4 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 18852-1-PBS in a different storage buffer formulation.