

For Research Use Only

# HHEX Polyclonal antibody

Catalog Number: 18882-1-AP



## Basic Information

Catalog Number: 18882-1-AP	GenBank Accession Number: BC014336	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 3087	
Source: Rabbit	Full Name: hematopoietically expressed homeobox	
Isotype: IgG	Calculated MW: 30 kDa	
Immunogen Catalog Number: AG13423		

## Applications

Tested Applications:  
IHC, ELISA

Species Specificity:  
human, mouse, rat

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data