

For Research Use Only

KLF1 Polyclonal antibody

Catalog Number: 18887-1-AP



Basic Information

Catalog Number:
18887-1-AP

Size:
150ul

Source:
Rabbit

Isotype:
IgG

Immunogen Catalog Number:
AG13467

GenBank Accession Number:
BC033580

GeneID (NCBI):
10661

Full Name:
Kruppel-like factor 1 (erythroid)

Calculated MW:
362 aa, 38 kDa

Purification Method:
Antigen affinity purification

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data