

For Research Use Only

PALM2 Polyclonal antibody

Catalog Number: 18889-1-AP



Basic Information

Catalog Number: 18889-1-AP	GenBank Accession Number: BC039306	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 700 µg/ml by Nanodrop and 553 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 114299	Recommended Dilutions: WB 1:5000-1:50000
Source: Rabbit	Full Name: paralemmin 2	
Isotype: IgG	Calculated MW: 1103aa,122 kDa; 377aa,42 kDa	
Immunogen Catalog Number: AG13471	Observed MW: 65-75 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, mouse brain tissue, rat brain tissue
Species Specificity: human, mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 18889-1-AP (PALM2 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.