

For Research Use Only

PALM2 Polyclonal antibody

Catalog Number:18889-1-AP



Basic Information

Catalog Number:

18889-1-AP

Size:

150ul , Concentration: 700 ug/ml by Nanodrop and 553 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13471

GenBank Accession Number:

BC039306

GeneID (NCBI):

114299

UNIPROT ID:

Q8IXS6

Full Name:

paralemmin 2

Calculated MW:

1103aa,122 kDa; 377aa,42 kDa

Observed MW:

65-75 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, mouse brain tissue, rat brain tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 18889-1-AP (PALM2 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.