

For Research Use Only

POU3F4-Specific Polyclonal antibody



Catalog Number: 19000-1-AP

Basic Information

Catalog Number: 19000-1-AP	GenBank Accession Number: NM_000307	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 747 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 5456	
Source: Rabbit	Full Name: POU class 3 homeobox 4	
Isotype: IgG	Calculated MW: 39 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

OTF9 also named as BRN4 and POU3F4, belongs to the POU transcription factor family and Class-3 subfamily. OTF9 is a transcription factor which exerts its primary action widely during early neural development and in a very limited set of neurons in the mature brain. Mutation of OTF9 will cause of X-linked deafness type 3 (DFN3). This antibody is specific to OTF9.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1%BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data