

For Research Use Only

H1FX Polyclonal antibody, PBS Only

Catalog Number:19321-1-PBS



Basic Information

Catalog Number:

19321-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5819

GenBank Accession Number:

BC000426

GeneID (NCBI):

8971

UNIPROT ID:

Q92522

Full Name:

H1 histone family, member X

Calculated MW:

22 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

H1FX (Histone H1x) belongs to the histone H1/H5 family which plays a crucial role in the organization and regulation of chromatin structure.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 19321-1-AP (H1FX antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 19321-1-PBS in a different storage buffer formulation.