

For Research Use Only

IL-10RB Polyclonal antibody

Catalog Number: 19387-1-AP



Basic Information

Catalog Number:

19387-1-AP

Size:

150ul , Concentration: 500 µg/ml by Nanodrop and 360 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7468

GenBank Accession Number:

BC001903

GeneID (NCBI):

3588

UNIPROT ID:

Q08334

Full Name:

interleukin 10 receptor, beta

Calculated MW:

37 kDa

Observed MW:

37-41 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Jurkat cells, HeLa cells, K-562 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 19387-1-AP (IL10RB antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.