

For Research Use Only

Aldehyde oxidase Polyclonal antibody, PBS Only

Catalog Number: 19495-1-PBS

Featured Product



Basic Information

Catalog Number:

19495-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001159

GeneID (NCBI):

316

UNIPROT ID:

Q06278

Full Name:

aldehyde oxidase 1

Calculated MW:

148 kDa

Observed MW:

148 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IF/ICC, IP, Indirect ELISA

Species Specificity:

human, mouse, rat, canine

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

