

For Research Use Only

# SP140 Polyclonal antibody

Catalog Number: 19506-1-AP



## Basic Information

<b>Catalog Number:</b> 19506-1-AP	<b>GenBank Accession Number:</b> NM_007237	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul	<b>GeneID (NCBI):</b> 11262	
<b>Source:</b> Rabbit	<b>Full Name:</b> SP140 nuclear body protein	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 98 kDa	

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human

## Background Information

SP140 is one component of the nuclear body, also known as nuclear domain 10, PML oncogenic domain, and KR body. SP140 may be involved in the pathogenesis of acute promyelocytic leukemia and viral infection. The antibody can recognize SP140 but can't recognize the isoform LYSP100-A.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com) W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data