

For Research Use Only

SP140 Polyclonal antibody

Catalog Number: 19506-1-AP



Basic Information

Catalog Number:

19506-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_007237

GeneID (NCBI):

11262

UNIPROT ID:

Q13342

Full Name:

SP140 nuclear body protein

Calculated MW:

98 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

SP140 is one component of the nuclear body, also known as nuclear domain 10, PML oncogenic domain, and KR body. SP140 may be involved in the pathogenesis of acute promyelocytic leukemia and viral infection. The antibody can recognize SP140 but can't recognize the isoform LYSP100-A.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

