

For Research Use Only

SPACA5 Polyclonal antibody

Catalog Number: 19525-1-AP



Basic Information

Catalog Number:

19525-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_205856

GeneID (NCBI):

389852

UNIPROT ID:

Q96QH8

Full Name:

sperm acrosome associated 5

Calculated MW:

18 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

SPACA5, also named as LYZL5, SPACA5A and SPACA5B, belongs to the glycosyl hydrolase 22 family. SPACA5 hydrolysis (1->4)-beta-linkages between N-acetylmuramic acid and N-acetyl-D-glucosamine residues in a peptidoglycan and between N-acetyl-D-glucosamine residues in chitodextrins.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

