

For Research Use Only

DNAI2 Polyclonal antibody

Catalog Number: 19541-1-AP



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
19541-1-AP	NM_023036	Antigen affinity purification
Size:	GeneID (NCBI):	
150ul	64446	
Source:	Full Name:	
Rabbit	dynein, axonemal, intermediate chain 2	
Isotype:	Calculated MW:	
IgG	69 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

DNAI2 belongs to the dynein intermediate chain family. It contains 5 WD repeats. DNAI2 is a part of the dynein complex of respiratory cilia. Mutation of DNAI2 may cause of primary ciliary dyskinesia type 9 (CILD9).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data