

For Research Use Only

# DNAI2 Polyclonal antibody

Catalog Number: 19541-1-AP



## Basic Information

Catalog Number:

19541-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_023036

GeneID (NCBI):

64446

UNIPROT ID:

Q9GZS0

Full Name:

dynein, axonemal, intermediate chain 2

Calculated MW:

69 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

## Background Information

DNAI2 belongs to the dynein intermediate chain family. It contains 5 WD repeats. DNAI2 is a part of the dynein complex of respiratory cilia. Mutation of DNAI2 may cause of primary ciliary dyskinesia type 9 (CILD9).

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

