

For Research Use Only

PLUNC Polyclonal antibody

Catalog Number: 19562-1-AP



Basic Information

| | | |
|--------------------------------------|---|--|
| Catalog Number: 19562-1-AP | GenBank Accession Number: NM_130852 | Purification Method: Antigen affinity purification |
| Size: 150ul | GeneID (NCBI): 51297 | |
| Source: Rabbit | Full Name: palate, lung and nasal epithelium associated | |
| Isotype: IgG | Calculated MW: 27 kDa | |

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

PLUNC, also named as LUNX, NASG and SPURT, belongs to the BPI/LBP/Plunc superfamily and Plunc family. PLUNC may be involved in the airway inflammatory response after exposure to irritants. It is associated with tumor progression. PLUNC plays a role in innate immune responses of the upper airways.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data