

For Research Use Only

PLUNC Polyclonal antibody

Catalog Number: 19562-1-AP



Basic Information

Catalog Number:

19562-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_130852

GeneID (NCBI):

51297

UNIPROT ID:

Q9NP55

Full Name:

palate, lung and nasal epithelium associated

Calculated MW:

27 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

PLUNC, also named as LUNX, NASG and SPURT, belongs to the BPI/LBP/Plunc superfamily and Plunc family. PLUNC may be involved in the airway inflammatory response after exposure to irritants. It is associated with tumor progression. PLUNC plays a role in innate immune responses of the upper airways.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

