

For Research Use Only

# ARHGEF10L Polyclonal antibody

Catalog Number: 19568-1-AP



## Basic Information

### Catalog Number:

19568-1-AP

### Size:

150ul , Concentration: 213 µg/ml by Bradford method using BSA as the standard;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG13746

### GenBank Accession Number:

BC003073

### GeneID (NCBI):

55160

### UNIPROT ID:

Q9HCE6

### Full Name:

Rho guanine nucleotide exchange factor (GEF) 10-like

### Calculated MW:

1279 aa, 140 kDa

### Purification Method:

Antigen affinity purification

## Applications

### Tested Applications:

ELISA

### Species Specificity:

human, mouse, rat

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

