

For Research Use Only

# CCT6B Polyclonal antibody

Catalog Number: 19668-1-AP



## Basic Information

<b>Catalog Number:</b> 19668-1-AP	<b>GenBank Accession Number:</b> NM_006584	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul	<b>GeneID (NCBI):</b> 10693	
<b>Source:</b> Rabbit	<b>Full Name:</b> chaperonin containing TCP1, subunit 6B (zeta 2)	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 58 kDa	

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human

## Background Information

CCT6B, also named as CCT-zeta-2, TCP-1-zeta-like, CCT-zeta-like, Testis-specific Tcp20 and Testis-specific protein TSA303, Belongs to the TCP-1 chaperonin family. CCT6B is a molecular chaperone which assists the folding of proteins upon ATP hydrolysis. It is known to play a role, in vitro, in the folding of actin and tubulin. The antibody has no cross reaction to CCT6A.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data