

For Research Use Only

CCT6B Polyclonal antibody

Catalog Number: 19668-1-AP



Basic Information

Catalog Number:

19668-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_006584

GeneID (NCBI):

10693

UNIPROT ID:

Q92526

Full Name:

chaperonin containing TCP1, subunit 6B (zeta 2)

Calculated MW:

58 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

CCT6B, also named as CCT-zeta-2, TCP-1-zeta-like, CCT-zeta-like, Testis-specific Tcp20 and Testis-specific protein TSA303, Belongs to the TCP-1 chaperonin family. CCT6B is a molecular chaperone which assists the folding of proteins upon ATP hydrolysis. It is known to play a role, in vitro, in the folding of actin and tubulin. The antibody has no cross reaction to CCT6A.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

