

For Research Use Only

DSG3-Specific Polyclonal antibody

Catalog Number: 19683-1-AP



Basic Information

Catalog Number: 19683-1-AP	GenBank Accession Number: NM_001944	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 220 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 1830	
Source: Rabbit	Full Name: desmoglein 3 (pemphigus vulgaris antigen)	
Isotype: IgG	Calculated MW: 108 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

DSG3, also named a sCDHF6 and PVA, is a component of intercellular desmosome junctions. DSG3 is involved in the interaction of plaque proteins and intermediate filaments mediating cell-cell adhesion. This antibody is specific to DSG3.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data