

For Research Use Only

ARSK-Specific Polyclonal antibody

Catalog Number: 19749-1-AP



Basic Information

Catalog Number:

19749-1-AP

Size:

150ul, Concentration: 900 ug/ml by Nanodrop and 580 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_198150

GeneID (NCBI):

153642

UNIPROT ID:

Q6UWY0

Full Name:

arylsulfatase family, member K

Calculated MW:

61 kDa

Observed MW:

55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells,

Background Information

ARSK, also named as TSULF and ASK, belongs to the sulfatase family. ARSK binds 1 calcium ion per subunit. The antibody is specific to ARSK.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 19749-1-AP (ARSK-Specific antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.