

For Research Use Only

ARSK-Specific Polyclonal antibody, PBS Only

Catalog Number: 19749-1-PBS



Basic Information

Catalog Number:
19749-1-PBS

Size:
100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:
Rabbit

Isotype:
IgG

GenBank Accession Number:
NM_198150

GeneID (NCBI):
153642

UNIPROT ID:
Q6UWY0

Full Name:
arylsulfatase family, member K

Calculated MW:
61 kDa

Observed MW:
55 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

ARSK, also named as TSULF and ASK, belongs to the sulfatase family. ARSK binds 1 calcium ion per subunit. The antibody is specific to ARSK.

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 19749-1-AP (ARSK-Specific antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 19749-1-PBS in a different storage buffer formulation.