

For Research Use Only

# OR10J5 Polyclonal antibody

Catalog Number: 19761-1-AP



## Basic Information

Catalog Number:

19761-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001004469

GeneID (NCBI):

127385

UNIPROT ID:

Q8NHC4

Full Name:

olfactory receptor, family 10,  
subfamily J, member 5

Calculated MW:

34 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

OR10J5, also named as OR1-28, belongs to the G-protein coupled receptor 1 family. OR10J5 is an odorant receptor.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

