

For Research Use Only

# VCAN Polyclonal antibody

Catalog Number: 19805-1-AP



## Basic Information

Catalog Number:

19805-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_004385

GeneID (NCBI):

1462

UNIPROT ID:

P13611

Full Name:

versican

Calculated MW:

373 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

VCAN, also named as PG-M, CSPG2 and GHAP, belongs to the aggrecan/versican proteoglycan family. VCAN may play a role in intercellular signaling and in connecting cells with the extracellular matrix. VCAN may take part in the regulation of cell motility, growth and differentiation. VCAN binds hyaluronic acid. The antibody recognizes the C-term of VCAN.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

