

For Research Use Only

VCAN Polyclonal antibody

Catalog Number: 19805-1-AP



Basic Information

Catalog Number: 19805-1-AP	GenBank Accession Number: NM_004385	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 1462	
Source: Rabbit	Full Name: versican	
Isotype: IgG	Calculated MW: 373 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

VCAN, also named as PG-M, CSPG2 and GHAP, belongs to the aggrecan/versican proteoglycan family. VCAN may play a role in intercellular signaling and in connecting cells with the extracellular matrix. VCAN may take part in the regulation of cell motility, growth and differentiation. VCAN binds hyaluronic acid. The antibody recognizes the C-term of VCAN.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data