

For Research Use Only

RGP1 Polyclonal antibody, PBS Only

Catalog Number:19908-1-PBS



Basic Information

Catalog Number:

19908-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13902

GenBank Accession Number:

BC001725

GeneID (NCBI):

9827

UNIPROT ID:

Q92546

Full Name:

RGP1 retrograde golgi transport homolog (S. cerevisiae)

Calculated MW:

391 aa, 42 kDa

Observed MW:

42 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 19908-1-AP (RGP1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 19908-1-PBS in a different storage buffer formulation.