

For Research Use Only

CD133-1,6,7 Polyclonal antibody

Catalog Number: 19944-1-AP **1 Publications**



Basic Information

Catalog Number: 19944-1-AP	GenBank Accession Number: NM_006017	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 8842	
Source: Rabbit	UNIPROT ID: O43490	
Isotype: IgG	Full Name: prominin 1	
	Calculated MW: 97 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human

Cited Species:
human, mouse, rat

Background Information

CD133, also known as PROM1 (prominin-1) or AC133, belongs to the prominin family. CD133 is a transmembrane glycoprotein with an NH₂-terminal extracellular domain, five transmembrane loops and a cytoplasmic tail. The expression of CD133 has been reported in hematopoietic stem cells, endothelial progenitor cells, neuronal and glial stem cells, suggesting the potential role of CD133 as a cell surface marker of adult stem cells. CD133 has also been reported as a marker of cancer stem cells in various human tumors. This antibody recognizes transcript variant 1,6,7 of CD133.

Notable Publications

Author	Pubmed ID	Journal	Application
Tao Li	41153008	Stem Cell Res Ther	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data