

For Research Use Only

CYP4F22-Specific Polyclonal antibody



Catalog Number: 20014-1-AP

Basic Information

Catalog Number:

20014-1-AP

Size:

150ul , Concentration: 173 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_173483

GeneID (NCBI):

126410

UNIPROT ID:

Q6NT55

Full Name:

cytochrome P450, family 4, subfamily F, polypeptide 22

Calculated MW:

62 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

CYP4F22, also named as FLJ39501 and LI3, belongs to the cytochrome P450 family. Cytochromes P450 are a group of heme-thiolate monooxygenases. Defects in CYP4F22 are the cause of ichthyosis lamellar type 3 (LI3). This antibody is specific to CYP4F22. It has no cross reaction to other CYP4F member.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

