

For Research Use Only

AKAP12 Polyclonal antibody

Catalog Number: 20016-1-AP



Basic Information

Catalog Number: 20016-1-AP	GenBank Accession Number: NM_005100	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 9590	
Source: Rabbit	Full Name: A kinase (PKA) anchor protein 12	
Isotype: IgG	Calculated MW: 191 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

AKAP12, also named as AKAP250, is an anchoring protein that mediates the subcellular compartmentation of protein kinase A (PKA) and protein kinase C (PKC). The antibody recognizes the C-term of AKAP12.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data