

For Research Use Only

KIAA0100 Polyclonal antibody

Catalog Number: 20019-1-AP



Basic Information

Catalog Number:

20019-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_014680

GeneID (NCBI):

9703

UNIPROT ID:

Q14667

Full Name:

KIAA0100

Calculated MW:

254 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

KIAA0100, also named as BCOX1 and MLAA-22, belongs to the UPF0378 family. The antibody recognizes all the isoform of KIAA0100.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

