

For Research Use Only

ANKRD26 Polyclonal antibody, PBS Only

Catalog Number: 20035-1-PBS



Basic Information

Catalog Number:

20035-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_014915

GeneID (NCBI):

22852

UNIPROT ID:

Q9UPS8

Full Name:

ankyrin repeat domain 26

Calculated MW:

196 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

ANKRD26 also named as KIAA1074. The antibody recognizes the C-term of ANKRD26.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human breast cancer slide using 20035-1-AP (ANKRD26 Antibody) at dilution of 1:50. This data was developed using the same antibody clone with 20035-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of HeLa cells using 20035-1-AP (ANKRD26 antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG. This data was developed using the same antibody clone with 20035-1-PBS in a different storage buffer formulation.