

For Research Use Only

HMGB3 Polyclonal antibody

Catalog Number: 20036-1-AP



Basic Information

Catalog Number:

20036-1-AP

Size:

150ul , Concentration: 253 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_014983

GeneID (NCBI):

22993

UNIPROT ID:

Q12766

Full Name:

HMG box domain containing 3

Calculated MW:

168 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

HMGB3 also named as KIAA0194 and SMF. The antibody recognizes near the N-term of HMGB3.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

