

For Research Use Only

PLXDC1 Polyclonal antibody

Catalog Number: 20051-1-AP



Basic Information

Catalog Number:

20051-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_020405

GeneID (NCBI):

57125

UNIPROT ID:

Q8IUK5

Full Name:

plexin domain containing 1

Calculated MW:

56 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

PLXDC1, also named as TEM3 and TEM7, belongs to the plexin family. PLXDC1 may interact with cortactin. PLXDC1(TEM7) is identified through serial analysis of gene expression on endothelial cells isolated from human normal and malignant colorectal tissues. The antibody recognizes isoform1 and isoform4 of PLXDC1.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

