

For Research Use Only

CDH24 Polyclonal antibody

Catalog Number: 20058-1-AP



Basic Information

Catalog Number: 20058-1-AP	GenBank Accession Number: NM_022478	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 64403	
Source: Rabbit	Full Name: cadherin-like 24	
Isotype: IgG	Calculated MW: 88 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

CDH24, also named as CDH11L, is a member of Cadherins which are calcium dependent cell adhesion proteins. Cadherine preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. Cadherin-24 mediates strong cell-cell adhesion. The antibody can recognize isoform1 and isoform2 of CDH24.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data