

For Research Use Only

CDH24 Polyclonal antibody

Catalog Number: 20058-1-AP



Basic Information

Catalog Number:

20058-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_022478

GeneID (NCBI):

64403

UNIPROT ID:

Q86UP0

Full Name:

cadherin-like 24

Calculated MW:

88 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

CDH24, also named as CDH11L, is a member of Cadherins which are calcium dependent cell adhesion proteins. Cadherine preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. Cadherin-24 mediates strong cell-cell adhesion. The antibody can recognize isoform1 and isoform2 of CDH24.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

